



## Marker Screening with VNR-145 Hybrid

Sr. No.	Hybrid	Marker	Polymorphic	GP Analysis
1	VNR-145	SSR-10	√	X
2	VNR-145	CASSR-22	√	√
3	VNR-145	SSR-116	X	X
4	VNR-145	CASSR-29	X	X
5	VNR-145	SSR-20	Not work	X

**Note:** Require parent samples for DNA fingerprinting to conformation & validation of Hybrid





**Report:** when I was screened some samples of VNR -145 with different markers than 2 type of polymorphism has been seen in SSR-10 & CaSSR-22. Some other markers also run but they were shown monomorphism in all samples, so am not include that marker in it.